



Technical Reproducibility Validation Report

STRsearch

Version c70179b

Autosomal STR

Verification date	2026-09-19
Repository commit	c70179b3b175adc82a7314409af06900b3861d61
Attestation level	Runs + Plausible output
Permalink	https://strhub.app/verified/strsearch
STRhub	https://strhub.app

SCOPE OF THIS REPORT

Independent automated verification of reproducible execution. It confirms the tool installs, runs end-to-end, and produces structurally valid output in a standardized environment.

This is **not** an analytical validation. Genotype accuracy, concordance, forensic suitability, and regulatory compliance are out of scope.

WHAT WAS VERIFIED

- Source available
- Installation successful
- End-to-end execution
- Output generated

NOT VERIFIED

- Genotype accuracy
- Concordance
- Forensic validity
- Regulatory compliance



1. Executive Summary

PURPOSE	Verify that the tool installs, runs end-to-end, and produces structurally valid output.
RESULT	5/5 gates passed. Runs + Plausible output
ATTESTATION LEVEL	Runs + Plausible output
SCOPE	Installation, execution, and output structure verification.
NOT EVALUATED	Genotype accuracy · Concordance · Forensic validity · Regulatory compliance

2. Reproducibility Metadata

Tool	STRsearch
Version	c70179b
Repository	https://github.com/AnJingwd/STRsearch
Commit (pinned)	c70179b3b175adc82a7314409af06900b3861d61
Container OS	ubuntu-22.04
Marker panel	Autosomal STR
Recipe	a recipe STRhub wrote by hand, not read off the repository: a note under the documented result, never the badge
Verified on	2026-09-19 18:20:00 UTC
CI run	https://github.com/Tfronta/strhub-verified/actions/runs/35460767865

Verified automatically, in a clean environment, on the tool's public source at the pinned commit. This certificate records what happened; it is not an endorsement by the tool's author.

3. Exact Run Command

This is the command that was executed, verbatim, in the CI environment.

```
# Environment: ubuntu-22.04 · STRsearch c70179b · commit c70179b
$ python3 /opt/STRsearch/pipeline.py \
--type paired \
--num_processors 4 \
--reads_threshold 30 \
--stutter_ratio 0.5 \
from_bam \
--working_path /data/out \
--sample sample \
--sex female \
--bam /data/in/input.bam \
--ref_bed /data/in/regions.bed \
--genotypes /data/out/genotypes.txt \
--multiple_alleles /data/out/multiple_alleles.txt \
--qc_matrix /data/out/qc_matrix.txt
```

4. Verification Gates

Available	The pinned public source exists	PASS
Installs	The environment builds from source	PASS
Runs	Executes end-to-end without crashing	PASS
Expected IO	Produces a non-empty file in the declared format	PASS
Plausible Output	Output structure matches expected genotype-bearing format	PASS
Execution log	strsearch.log-external.txt	

5. Out of Scope

This report does not evaluate any of the following:

- Genotype correctness or accuracy
- Concordance against known truth sets
- Sensitivity, specificity, or stutter performance
- Allele calling accuracy or forensic casework suitability
- Regulatory compliance or ISO accreditation
- Multi-laboratory or multi-dataset reproducibility

6. Output Content Evidence

Output file	multiple_alleles.txt
Format	TSV
Sequence records	388
Distinct STR loci	24
Total reads	2,401
Max depth (single locus)	51

7. Verification Data

The repository ships example data (35 file(s)); this run did not use it. STRhub ran the verification using the public reference dataset listed below. That dataset is a slice around 24 forensic STR loci, not a whole genome: it carries reads only at the loci listed below. A small test file in the repository lets a new user run the tool on their first day and see it working before trusting it with their own data, and it lets a verification run against the author's own sample as well as this one. Publishing the output that file should produce helps just as much: it shows what the results are meant to look like, which is what a reader needs to tell a correct run from one that merely finished.

Dataset	Illumina BAM (hg38), NA12878 (autosomal, female)
Source	https://ftp-trace.ncbi.nlm.nih.gov/ReferenceSamples/giab/dat...
License	Open access (GIAB / NIST). Research and educational use. STRhub is not the data provider.
Ref. genome	GRCh38 / hg38
Loci tested	CSF1PO · D10S1248 · D12S391 · D13S317 · D16S539 · D18S51 · D19S433 · D1S1656 · D21S11 · D22S1045 · D2S1338 · D2S441 · D3S1358 · D5S818 · D6S1043 · D7S820 · D8S1179 · FGA · PentaD · PentaE · SE33 · TH01 · TPOX · vWA
Regions BED	Provided by STRhub

8. Verification Matrix

PASS	Reference dataset	Illumina BAM (hg38), NA12878 (autosomal, female)
	README check (5/5)	Advisory: all items present
	Install / environment setup	Present
	Run command	Present
	Expected input format	Present
	Produced output format	Present
	Dependencies / versions	Present

9. What This Run Needed Beyond the Repository

The result above describes a run configured as follows. Anyone repeating it needs the same things.

- Test data: no sample from the repository was used, so a public reference sample stood in.
- A container environment, supplied with the submission.

10. What STRhub's recipe does that the repository's instructions do not

This run used a recipe STRhub wrote. Each item below is a departure from the README that a first-time user would have to discover for themselves. The badge rests on the run of the repository's own instructions, not on this one.

- Builds the environment on ubuntu 22.04 at the version floors the README declares (bwa, samtools, bamToFastq, seqtk). — instead of: The repository's own Dockerfile: ubuntu 16.04 with exact pins and Miniconda 4.3.31. The pinned stack no longer resolves.
- Fills in conf.py with the paths of the installed tools. — instead of: conf.py as shipped, with 'YOUR PATH/...' placeholders the user is to replace.
- Leaves usearch out. — instead of: usearch, which conf.py names. Licensed and unobtainable; only --assemble_pairs reaches it, and this run does not use it.
- Supplies an hg38 regions configuration in STRsearch's eleven-column layout, with flanking sequences, written by STRhub. — instead of: The repository's example, which is hg19 and covers five markers. No hg38 configuration ships with the repository.

11. Limitations

- Single reference dataset per input type
- Single containerized environment (Docker / ubuntu-22.04)
- No truth-set comparison or ground-truth genotypes
- No accuracy or concordance assessment
- No forensic validation of results
- Short-read limitations apply (very long STR alleles may not span reads)

12. Scope and Disclaimers

Executed end-to-end in the stated environment with output in the expected format. Concerns reproducible execution only; no claim of accuracy, casework fitness, or regulatory validation.

Each result is a dated snapshot, verified on the tool's public repository at a pinned commit. STRhub does not store tool source code.

13. Conclusion

◆ Runs end-to-end

STRsearch c70179b installs and executes without error in a clean ubuntu-22.04 environment at the pinned commit.

◆ Produces plausible output

The tool generated a structurally valid TSV output file with genotype calls across 24 target forensic STR loci, with a maximum read depth of 51 at vWA.

◆ Verdict: Runs

The tool installed and its run produced its documented output.

◆ Repository test data not used

The repository ships example data (35 file(s)), which this run did not use: the gates above ran on STRhub's reference dataset. STRhub Verified supplied Illumina BAM (hg38), NA12878 (autosomal, female) as the reference input for this verification run. Users replicating this result should use the same or equivalent publicly available reference material.

◆ Reproducibility statement

The exact command reported in Section 3, executed at the pinned commit in the described environment, is sufficient to reproduce this verification result independently.



LOCUS	READ DEPTH		LOCUS	READ DEPTH	
vWA	241	████████████████████	CSF1PO	92	██████
D21S11	187	██████████████	D6S1043	92	██████
FGA	164	██████████	D2S441	88	████
D19S433	131	██████	D1S1656	81	████
SE33	125	██████	PentaE	78	████
D18S51	123	██████	D13S317	66	███
D16S539	105	████	D10S1248	65	███
D3S1358	99	████	D2S1338	62	███
PentaD	99	████	D8S1179	62	███
D7S820	96	████	TPOX	55	██
D22S1045	95	████	D12S391	55	██
D5S818	94	████	TH01	46	██